



DNAAF3 Antibody, Biotin conjugated

Product Code	CSB-PA850845LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N9W5
Immunogen	Recombinant Human Dynein assembly factor 3, axonemal protein (1-209AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dynein assembly factor 3, axonemal, DNAAF3, C19orf51
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DNAAF3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.