



SLC16A10 Antibody, HRP conjugated

Product Code	CSB-PA851563HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TF71
Immunogen	Recombinant Human Monocarboxylate transporter 10 protein (159-201AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Monocarboxylate transporter 10 (MCT 10) (Aromatic amino acid transporter 1) (Solute carrier family 16 member 10) (T-type amino acid transporter 1), SLC16A10, MCT10 TAT1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC16A10