



Pdgfd Antibody, HRP conjugated

Product Code	CSB-PA852889LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q925I7
Immunogen	Recombinant Mouse Platelet-derived growth factor D protein (24-370AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Platelet-derived growth factor D (PDGF-D) (Spinal cord-derived growth factor B) (SCDGF-B) [Cleaved into: Platelet-derived growth factor D, latent form (PDGFD latent form); Platelet-derived growth factor D, receptor-binding form (PDGFD receptor-binding form)], Pdgfd, Scdgfb
Immunogen Species	Mus musculus (Mouse)
Research Area	Cardiovascular
Target Names	Pdgfd
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.