





EZH1 Antibody

Product Code	CSB-PA852913DA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92800
Immunogen	Recombinant Human Histone-lysine N-methyltransferase EZH1 protein (160-280AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Histone-lysine N-methyltransferase EZH1 (EC 2.1.1.43) (ENX-2) (Enhancer of zeste homolog 1), EZH1, KIAA0388
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	EZH1

Image



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA852913DA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA852913DA01HU at dilution of 1:100