



NEO1 Antibody, FITC conjugated

Product Code	CSB-PA852917LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92859
Immunogen	Recombinant Human Neogenin protein (414-682AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neogenin (Immunoglobulin superfamily DCC subclass member 2), NEO1, IGDC2 NGN
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NEO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.