



MR1 Antibody, Biotin conjugated

Product Code	CSB-PA853268LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q95460
Immunogen	Recombinant Human Major histocompatibility complex class I-related gene protein (150-253AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major histocompatibility complex class I-related gene protein (MHC class I-related gene protein) (Class I histocompatibility antigen-like protein), MR1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	MR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.