




TMC2 Antibody

Product Code	CSB-PA855056LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDI7
Immunogen	Recombinant Human Transmembrane channel-like protein 2 protein (6-206AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane channel-like protein 2 (Transmembrane cochlear-expressed protein 2), TMC2, C20orf145
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TMC2

Image



Immunofluorescent analysis of HeLa cells using CSB-PA855056LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.