



TMC2 Antibody, FITC conjugated

Product Code	CSB-PA855056LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDI7
Immunogen	Recombinant Human Transmembrane channel-like protein 2 protein (6-206AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane channel-like protein 2 (Transmembrane cochlear-expressed protein 2), TMC2, C20orf145
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TMC2