



EAF2 Antibody, FITC conjugated

Product Code	CSB-PA856910LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96CJ1
Immunogen	Recombinant Human ELL-associated factor 2 protein (6-100AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ELL-associated factor 2 (Testosterone-regulated apoptosis inducer and tumor suppressor protein), EAF2, TRAITS
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	EAF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.