



PEBP4 Antibody, Biotin conjugated

Product Code	CSB-PA857025LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96S96
Immunogen	Recombinant Human Phosphatidylethanolamine-binding protein 4 protein (23-227AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylethanolamine-binding protein 4 (PEBP-4) (hPEBP4) (Protein cousin-of-RKIP 1), PEBP4, CORK1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PEBP4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.