



ARHGDIG Antibody, HRP conjugated

Product Code	CSB-PA858727LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99819
Immunogen	Recombinant Human Rho GDP-dissociation inhibitor 3 protein (1-225AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho GDP-dissociation inhibitor 3 (Rho GDI 3) (Rho-GDI gamma), ARHGDIG
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARHGDIG