



PIF1 Antibody, Biotin conjugated

Product Code	CSB-PA863987LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H611
Immunogen	Recombinant Human ATP-dependent DNA helicase PIF1 protein (125-360AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-dependent DNA helicase PIF1 (EC 3.6.4.12) (DNA repair and recombination helicase PIF1) (PIF1/RRM3 DNA helicase-like protein), PIF1, C15orf20
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PIF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.