





DPH7 Antibody

Product Code	CSB-PA866281LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BTV6
Immunogen	Recombinant Human Diphthine methyltransferase protein (210-325AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Diphthine methyltransferase (EC 3.1.1.97) (Diphthamide biosynthesis protein 7) (WD repeat-containing protein 85), DPH7, C9orf112 WDR85
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DPH7


Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA866281LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA866281LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA866281LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)