



# DPH7 Antibody, Biotin conjugated

Product Code	CSB-PA866281LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BTV6
Immunogen	Recombinant Human Diphthine methyltransferase protein (210-325AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Diphthine methyltransferase (EC 3.1.1.97) (Diphthamide biosynthesis protein 7) (WD repeat-containing protein 85), DPH7, C9orf112 WDR85
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DPH7