



DPH7 Antibody, Biotin conjugated

Product Code	CSB-PA866281LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BTV6
Immunogen	Recombinant Human Diphthine methyltransferase protein (210-325AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Diphthine methyltransferase (EC 3.1.1.97) (Diphthamide biosynthesis protein 7) (WD repeat-containing protein 85), DPH7, C9orf112 WDR85
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DPH7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.