



SCPEP1 Antibody, HRP conjugated

Product Code	CSB-PA867189LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HB40
Immunogen	Recombinant Human Retinoid-inducible serine carboxypeptidase protein (201-452AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Retinoid-inducible serine carboxypeptidase (EC 3.4.16.-) (Serine carboxypeptidase 1), SCPEP1, RISC SCP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SCPEP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.