



SPATA16 Antibody, FITC conjugated

Product Code	CSB-PA871597LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXB7
Immunogen	Recombinant Human Spermatogenesis-associated protein 16 protein (238-475AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Spermatogenesis-associated protein 16 (Testis development protein NYD-SP12), SPATA16
Immunogen Species	Homo sapiens (Human)
Research Area	Developmental Biology
Target Names	SPATA16