



IFT46 Antibody, Biotin conjugated

Product Code	CSB-PA873621LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQC8
Immunogen	Recombinant Human Intraflagellar transport protein 46 homolog protein (1-304AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Intraflagellar transport protein 46 homolog, IFT46, C11orf2 C11orf60
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	IFT46
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.