



# EXD2 Antibody, Biotin conjugated

Product Code	CSB-PA873662LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NVH0
Immunogen	Recombinant Human Exonuclease 3'-5' domain-containing protein 2 protein (137-245AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Exonuclease 3'-5' domain-containing protein 2 (EC 3.1.11.1) (Exonuclease 3'-5' domain-like-containing protein 2), EXD2, C14orf114 EXDL2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	EXD2