



PLAC8 Antibody, HRP conjugated

Product Code	CSB-PA873705LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZF1
Immunogen	Recombinant Human Placenta-specific gene 8 protein (1-110AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Placenta-specific gene 8 protein (Protein C15), PLAC8
Immunogen Species	Homo sapiens (Human)
Research Area	Developmental Biology
Target Names	PLAC8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.