



TDRD1 Antibody, Biotin conjugated

Product Code	CSB-PA874834LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXT4
Immunogen	Recombinant Human Tudor domain-containing protein 1 protein (12-314AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tudor domain-containing protein 1 (Cancer/testis antigen 41.1) (CT41.1), TDRD1
Immunogen Species	Homo sapiens (Human)
Research Area	Developmental Biology
Target Names	TDRD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.