



SMOC1 Antibody, HRP conjugated

Product Code	CSB-PA875673LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4F8
Immunogen	Recombinant Human SPARC-related modular calcium-binding protein 1 protein (277-382AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SPARC-related modular calcium-binding protein 1 (Secreted modular calcium-binding protein 1) (SMOC-1), SMOC1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SMOC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.