



GMIP Antibody, Biotin conjugated

Product Code	CSB-PA878952LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P107
Immunogen	Recombinant Human GEM-interacting protein (719-858AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GEM-interacting protein (GMIP), GMIP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GMIP