



DNAH1 Antibody, Biotin conjugated

Product Code	CSB-PA878961LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P2D7
Immunogen	Recombinant Human Dynein heavy chain 1, axonemal protein (2225-2407AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dynein heavy chain 1, axonemal (Axonemal beta dynein heavy chain 1) (Ciliary dynein heavy chain 1) (Heat shock regulated protein 1) (HSRF-1) (hDHC7), DNAH1, DHC7 DNAHC1 KIAA1410
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DNAH1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.