




BCL11B Antibody

Product Code	CSB-PA880159HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9C0K0
Immunogen	Recombinant Human B-cell lymphoma/leukemia 11B protein (500-550AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:200-1:500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	B-cell lymphoma/leukemia 11B (BCL-11B) (B-cell CLL/lymphoma 11B) (COUP-TF-interacting protein 2) (Radiation-induced tumor suppressor gene 1 protein) (hRit1), BCL11B, CTIP2 RIT1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	BCL11B

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA880159HA01HU at dilution of 1:300

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.