




GBA2 Antibody

Product Code	CSB-PA881025LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HCG7
Immunogen	Recombinant Human Non-lysosomal glucosylceramidase protein (656-867AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Non-lysosomal glucosylceramidase (NLGase) (EC 3.2.1.45) (Beta-glucocerebrosidase 2) (Beta-glucosidase 2) (Glucosylceramidase 2), GBA2, KIAA1605
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GBA2

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA881025LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.