



GBA2 Antibody, Biotin conjugated

Product Code	CSB-PA881025LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HCG7
Immunogen	Recombinant Human Non-lysosomal glucosylceramidase protein (656-867AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Non-lysosomal glucosylceramidase (NLGase) (EC 3.2.1.45) (Beta-glucocerebrosidase 2) (Beta-glucosidase 2) (Glucosylceramidase 2), GBA2, KIAA1605
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GBA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.