



FLRT3 Antibody, HRP conjugated

Product Code	CSB-PA882162LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZU0
Immunogen	Recombinant Human Leucine-rich repeat transmembrane protein FLRT3 protein (209-264AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich repeat transmembrane protein FLRT3 (Fibronectin-like domain-containing leucine-rich transmembrane protein 3), FLRT3, KIAA1469
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FLRT3