



SLC25A10 Antibody, FITC conjugated

Product Code	CSB-PA883371LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBX3
Immunogen	Recombinant Human Mitochondrial dicarboxylate carrier protein (124-162AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial dicarboxylate carrier (Solute carrier family 25 member 10), SLC25A10, DIC
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC25A10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.