



BAZ2A Antibody, FITC conjugated

Product Code	CSB-PA883428HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UIF9
Immunogen	Recombinant Human Bromodomain adjacent to zinc finger domain protein 2A protein (97-266AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bromodomain adjacent to zinc finger domain protein 2A (Transcription termination factor I-interacting protein 5) (TTF-I-interacting protein 5) (Tip5) (hWALp3), BAZ2A, KIAA0314 TIP5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	BAZ2A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.