



BAZ2A Antibody, FITC conjugated

Product Code	CSB-PA883428HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UIF9
Immunogen	Recombinant Human Bromodomain adjacent to zinc finger domain protein 2A protein (97-266AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bromodomain adjacent to zinc finger domain protein 2A (Transcription termination factor I-interacting protein 5) (TTF-I-interacting protein 5) (Tip5) (hWALp3), BAZ2A, KIAA0314 TIP5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	BAZ2A