




CYREN Antibody

Product Code	CSB-PA883613HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BWK5
Immunogen	Recombinant Human Cell cycle regulator of non-homologous end joining protein (1-157AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell cycle regulator of non-homologous end joining, Cell cycle regulator of NHEJ, Modulator of retrovirus infection homolog, CYREN, MRI, C7orf49
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CYREN

Image



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA883613HA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.