



CYREN Antibody, Biotin conjugated

Product Code	CSB-PA883613HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BWK5
Immunogen	Recombinant Human Cell cycle regulator of non-homologous end joining protein (1-157AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell cycle regulator of non-homologous end joining, Cell cycle regulator of NHEJ, Modulator of retrovirus infection homolog, CYREN, MRI, C7orf49
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CYREN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.