



CYREN Antibody, Biotin conjugated

Product Code	CSB-PA883613HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BWK5
Immunogen	Recombinant Human Cell cycle regulator of non-homologous end joining protein (1-157AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell cycle regulator of non-homologous end joining, Cell cycle regulator of NHEJ, Modulator of retrovirus infection homolog, CYREN, MRI, C7orf49
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CYREN