




EHD1 Antibody

Product Code	CSB-PA884470LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4M9
Immunogen	Recombinant Human EH domain-containing protein 1 protein (401-498AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	EH domain-containing protein 1 (PAST homolog 1) (hPAST1) (Testilin), EHD1, PAST PAST1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	EHD1

Image



Western blot


All lanes: EHD1 antibody at 12µg/ml + MDA-MB-231 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 61 kDa

Observed band size: 61 kDa



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA884470LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.