





# EHD1 Antibody

<b>Product Code</b>	CSB-PA884470LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9H4M9
<b>Immunogen</b>	Recombinant Human EH domain-containing protein 1 protein (401-498AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	EH domain-containing protein 1 (PAST homolog 1) (hPAST1) (Testilin), EHD1, PAST PAST1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cardiovascular
<b>Target Names</b>	EHD1

## Image



Western blot  
All lanes: EHD1 antibody at 12µg/ml + MDA-MB-231 whole cell lysate  
Secondary  
Goat polyclonal to rabbit IgG at 1/10000 dilution  
Predicted band size: 61 kDa  
Observed band size: 61 kDa



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA884470LA01HU at dilution of 1:100