



EHD1 Antibody, Biotin conjugated

Product Code	CSB-PA884470LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4M9
Immunogen	Recombinant Human EH domain-containing protein 1 protein (401-498AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	EH domain-containing protein 1 (PAST homolog 1) (hPAST1) (Testilin), EHD1, PAST PAST1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	EHD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.