



EHD1 Antibody, Biotin conjugated

Product Code	CSB-PA884470LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4M9
Immunogen	Recombinant Human EH domain-containing protein 1 protein (401-498AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	EH domain-containing protein 1 (PAST homolog 1) (hPAST1) (Testilin), EHD1, PAST PAST1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	EHD1