




EBF2 Antibody

Product Code	CSB-PA884503LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HAK2
Immunogen	Recombinant Human Transcription factor COE2 protein (53-254AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor COE2 (Early B-cell factor 2) (EBF-2), EBF2, COE2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	EBF2

Image



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA884503LA01HU at dilution of 1:100