



OVOL2 Antibody, FITC conjugated

Product Code	CSB-PA887110HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRP0
Immunogen	Recombinant Human Transcription factor Ovo-like 2 protein (19-114AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor Ovo-like 2 (hOvo2) (Zinc finger protein 339), OVOL2, ZNF339
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	OVOL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.