



TBC1D24 Antibody, Biotin conjugated

Product Code	CSB-PA890752LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULP9
Immunogen	Recombinant Human TBC1 domain family member 24 protein (1-169AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TBC1 domain family member 24, TBC1D24, KIAA1171
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TBC1D24
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.