



# FBXL21 Antibody, Biotin conjugated

Product Code	CSB-PA891953LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKT6
Immunogen	Recombinant Human F-box/LRR-repeat protein 21 protein (78-137AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box/LRR-repeat protein 21 (F-box and leucine-rich repeat protein 21) (F-box and leucine-rich repeat protein 3B) (F-box/LRR-repeat protein 3B), FBXL21, FBL21 FBL3 FBXL3B FBXL3P
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FBXL21