



PPM1H Antibody, Biotin conjugated

Product Code	CSB-PA892345LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULR3
Immunogen	Recombinant Human Protein phosphatase 1H protein (73-289AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein phosphatase 1H (EC 3.1.3.16), PPM1H, ARHCL1 KIAA1157 URCC2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PPM1H