



PPM1H Antibody, Biotin conjugated

Product Code	CSB-PA892345LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULR3
Immunogen	Recombinant Human Protein phosphatase 1H protein (73-289AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein phosphatase 1H (EC 3.1.3.16), PPM1H, ARHCL1 KIAA1157 URCC2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PPM1H
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.