





DCAF1 Antibody

Product Code	CSB-PA896722LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y4B6
Immunogen	Recombinant Human DDB1- and CUL4-associated factor 1 protein (647-745AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DDB1- and CUL4-associated factor 1 (HIV-1 Vpr-binding protein) (VprBP) (Serine/threonine-protein kinase VPRBP) (EC 2.7.11.1) (Vpr-interacting protein), DCAF1, KIAA0800 RIP VPRBP
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	DCAF1


Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA896722LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA896722LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA896722LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.