



ICE1 Antibody, Biotin conjugated

Product Code	CSB-PA897466LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2F5
Immunogen	Recombinant Human Little elongation complex subunit 1 protein (1213-1505AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Little elongation complex subunit 1 (Interactor of little elongator complex ELL subunit 1), ICE1, KIAA0947
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ICE1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.