



MAST1 Antibody, FITC conjugated

Product Code	CSB-PA897529LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2H9
Immunogen	Recombinant Human Microtubule-associated serine/threonine-protein kinase 1 protein (687-879AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Microtubule-associated serine/threonine-protein kinase 1 (EC 2.7.11.1) (Syntrophin-associated serine/threonine-protein kinase), MAST1, KIAA0973 SAST
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MAST1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.