




Recombinant Human Ras association domain-containing protein 6 (RASSF6)

Product Code	CSB-EP754570HU
Relevance	Involved in the induction of apoptosis, through both caspase-dependent and caspase-independent pathways. May act as a Ras effector protein. May suppress the serum-induced basal levels of NF-kappa-B
Abbreviation	Recombinant Human RASSF6 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZTQ3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MTMMAHQYPSWIFINEKTFITREQLNSLLKTYNIFYENQKNLHILYGETEDGKLI VEGMLDIFWGVKRPIQLKIQDEKPFSSFTSMKSSDVFSSKGMTRWGEFDDLY RISELDRTQIPMSEKRNSQEDYLSYHSNTLKPFAKDEPDSPVLYRTMSEALV RKRMPKPLMMDRKERQKNRASINGHFYNHEIENSPQDFALHIIFATGEQRRLKK TDIPLLQRLLQGPSEKNARIFLMDKDAEEISSDVAQYINFHFSLLLESILQRLNEEE KREIQRIVTKFNKEKAILKCLQNKLVIKTETTV
Research Area	others
Source	E.coli
Target Names	RASSF6
Protein Names	Recommended name: Ras association domain-containing protein 6
Expression Region	1-303aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	62.8kDa
Protein Length	Full Length of Isoform 4
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.