




Recombinant Human Tripartite motif-containing protein 69 (TRIM69)

Product Code	CSB-EP774824HU
Relevance	May have E3 ubiquitin-protein ligase activity. May play a role in apoptosis.
Abbreviation	Recombinant Human TRIM69 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86WT6
Alias	RFP-like domain-containing protein trimless RING finger protein 36 Tripartite motif-containing protein 69
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MEEELAIQQGQLETTLKEQLTLRNMQKEAIAAHKENKLHLQQHVSMEFLKLHQ FLHSKEKDILTELREEGKALNEEMELNLSQLQEQLLAKDMLVSIQAKTEQQNS FDFLKDITLLHSLEQGMKVLATRELISRKLNLGQYKGPQYQMVVWREMQDTLCP GLSPLTLDPKTAHPNLVLSKSQTSVWHGDIKKIMPDDPERFDSSVAVLGSRGF TSGKQWYWEVEVAKKTKWTVGVVRESIIRKGCPLTPEQGFWLLRLRNQTDLK ALDLPSFSLTLTNNLDKVGIYLDYEGGQLSFYNAKTMTHIYTFSTNFMEKLYPY FCPCCLNDGGENKEPLHILHPQ
Research Area	Cell Biology
Source	E.coli
Target Names	TRIM69
Protein Names	Recommended name: Tripartite motif-containing protein 69 Alternative name(s): RFP-like domain-containing protein trimless RING finger protein 36
Expression Region	1-341aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	66.3kDa
Protein Length	Full Length of Isoform 2
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.