




# Recombinant Mouse Trifunctional enzyme subunit alpha, mitochondrial (Hadha), partial

<b>Product Code</b>	CSB-BP804834MO
<b>Abbreviation</b>	Recombinant Mouse Hadha protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8BMS1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥ 85% as determined by SDS-PAGE.
<b>Sequence</b>	LCKKNKFGAPQKNVQQLAILGAGLMGAGIAQVSVDKGLKTLKDTTDTVTGLGRG QQQVFKGLNDKVKKKALTSFERDSIFSNLIGQLDYKGFEEKADMVIEAVFEDLGV KHKVLKEVESVTPEHCIFASNTSALPINQIAAVSKRPEKVGIMHYFSPVDKMQLL EIITTDKTSKDTTASAVAVGLRQQKVIVVKDGPFGFYTTTRCLAPMMSEVMRILQE GVDPKKLDALTTGFGFPVGAATLADEVGVDVAQHVAEDLGKAFGERFGGGSV ELLKQMVSKGFLGRKSGKGFYIYQEGSKNKSLSNEMDNILANLRLPAKPEVSS DEDVQYRVITRFVNEAVLCLQEGILATPAEGDIGAVFGLGFPPCLGGPFRFVDL YGAQKVVDRLRKYESAYGTQFTPCQLLLDHANNSSKKFYQ
<b>Research Area</b>	Cancer
<b>Source</b>	Baculovirus
<b>Target Names</b>	Hadha
<b>Protein Names</b>	Recommended name: Trifunctional enzyme subunit alpha, mitochondrial Alternative name(s): TP-alpha Including the following 2 domains: Long-chain enoyl-CoA hydratase EC= 4.2.1.17 Long chain 3-hydroxyacyl-CoA dehydrogenase EC= 1.1.1
<b>Expression Region</b>	348-763aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged and C-terminal Myc-tagged
<b>Mol. Weight</b>	50.1 kDa
<b>Protein Length</b>	Partial
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.